


1873 NECK STERILIZER SIDEFLEXING CHAIN WITH BASE ROLLER CHAIN

TARGET INDUSTRIES:

BEVERAGE: NECK/CAP STERILIZING HOT FILL

The product is specially designed to tilt hot-filled bottles on their side to allow hot liquid to sterilize the inside of the caps and then put them back upright.





ADVANTAGES

- The special design allows for gentle handling of the bottles and avoids bottle distortion and extra pressure on the bottles.
- The design of the rubber on top of the plate gives the bottle a good and stable support
- Flexibility of the chain for compact design if required
- Materials that can handle high temperatures
- Plate in low friction material
- Flexible base chain, can twist 90° in 1500mm.
- Metal base chain provides increased strength



- Longer life
- Lower friction
- Good chemical resistance
- High abrasion resistance

- Approved for direct food contact
- Less dust than with acetal chainstiles and avoids bottle distortion and extra pressure on the bottles.



**SYSTEM
PLAST™**

System Plast S.R.L.

Via Guareschi, nr. 2, 24060 Telgate (BG)
 Customer Service: +39 035. 835 1301
 Fax: +39 053. 835 1399
 E-mail: info@systemplast.com

www.RegalPTS.com

APPLICATION CONSIDERATIONS

The proper selection and application of power transmission products and components, including the related area of product safety, is the responsibility of the customer. Operating and performance requirements and potential associated issues will vary appreciably depending upon the use and application of such products and components. The scope of the technical and application information included in this publication is necessarily limited. Unusual operating environments and conditions, lubrication requirements, loading supports, and other factors can materially affect the application and operating results of the products and components and the customer should carefully review its requirements. Any technical advice or review furnished by Regal Beloit America, Inc. and its affiliates with respect to the use of products and components is given in good faith and without charge, and Regal assumes no obligation or liability for the advice given, or results obtained, all such advice and review being given and accepted at customer's risk. For a copy of our Standard Terms and Conditions of Sale, Disclaimers of Warranty, Limitation of Liability and Remedy, please contact Customer Service at 1-800-626-2120. These terms and conditions of sale, disclaimers and limitations of liability apply to any person who may buy, acquire or use a Regal Beloit America Inc. product referred to herein, including any person who buys from a licensed distributor of these branded products.

System Plast and Regal are trademarks of Regal Beloit Corporation or one of its affiliated companies.
 ©2016, 2016 Regal Beloit Corporation, All Rights Reserved. MCSSP16013E • Form 9974E • Printed in USA

Product name	Description
Sideflexing Chain	NGE1873SS-TAB-K1161VG85H8.5

*1873 Neck Sterilizer Sideflexing Chain with Base Roller Chain is not available in Germany

